

SpeedySwab Rapid COVID-19 + FLU A&B

Up To 15 Month Shelf Life

AVAILABLE IN 2 TEST-PACK



Convenient, Reliable At-Home Testing

The ability to quickly and accurately differentiate between COVID-19 and flu symptoms is more crucial than ever. This innovative Combo COVID-19 and Flu OTC Test empowers individuals with a convenient, reliable solution for at-home testing.



C = Control Line
 A = Influenza A Line
 B = Influenza B Line
 S = COVID-19 (SARS-CoV 2)

Who is it for?

- Individuals with signs and symptoms of respiratory infection consistent with COVID-19 within the first five (5) days of symptom onset
- Individuals aged 14 years or older will self-collect anterior nasal swab specimen
- Individuals age two (2) years or older will have adult-collect anterior nasal swab specimen

Quick and Easy Workflow

01 GENTLY SWAB AND COLLECT

02 INSERT SWAB AND STIR

03 APPLY SAMPLE ON TEST DEVICE

04 WAIT

READ RESULTS

SpeedySwab Rapid COVID-19 + FLU A&B

Up To 15 Month Shelf Life

AVAILABLE IN 2 TEST-PACK



<p>Intended Use</p>	<p>The SpeedySwab™ Rapid COVID-19 + Flu A&B Antigen Self-Test is a lateral flow immunoassay intended for the qualitative detection and differentiation of SARS-CoV-2, influenza A, and influenza B protein antigens. This test is authorized for non-prescription home use with self-collected anterior nasal swab specimens from individuals aged 14 years or older, or with adult-collected anterior nasal swab specimens from individuals two (2) years or older. This test is only authorized for individuals with signs and symptoms of respiratory infection consistent with COVID-19 within the first five (5) days of symptom onset when tested at least twice over three days with at least 48 hours between tests.</p>
<p>Sample Type</p>	<p>Nasal Swab</p>
<p>Result Time</p>	<p>15 Minutes</p>
<p>Clinical Sensitivity</p>	<p>COVID - 92.6%, FLU A - 82.9%, FLU B - 90%</p>
<p>Test Components</p>	<p>Sterile Swab x2, Test Tube with Buffer Solution x2, Test Device x2, Tube Holder x1, Quick Reference Instructions</p>
<p>EUA Number</p>	<p>EUA240014</p>