


SpeedySwab Rapid COVID-19 + FLU A&B

Up To 18 Month Shelf Life

AVAILABLE IN 25 TEST-PACK



Rapid Results, Batch Testing

The ability to quickly and accurately differentiate between COVID-19 and flu symptoms is more crucial than ever. This innovative Combo COVID-19 and Flu POC Test empowers providers with a convenient and reliable solution for in-clinic testing.




- C = Control Line
- A = Influenza A Line
- B = Influenza B Line
- S = COVID-19 (SARS-CoV 2)

Who is it for?

Individuals with signs and symptoms of respiratory infection consistent with COVID-19 within the first five (5) days of symptom onset


Quick and Easy Workflow



SpeedySwab Rapid COVID-19 + FLU A&B

Up To 18 Month Shelf Life

AVAILABLE IN 25 TEST-PACK



Intended Use	The SpeedySwab™ Rapid COVID-19 + Flu A&B Antigen Test is a lateral flow immunochromatographic assay intended for in-vitro rapid, simultaneous qualitative detection and differentiation of influenza A and B nucleoprotein antigens and SARS-CoV-2 nucleocapsid antigen directly from anterior nasal swab specimens of individuals with signs and symptoms of respiratory infection consistent with COVID-19 within the first five (5) days of symptom onset when tested at least twice over three days with at least 48 hours between tests.
Sample Type	Nasal Swab
Result Time	15 Minutes
Clinical Sensitivity	COVID - 92.6%, FLU A - 82.9%, FLU B - 90%
Test Components	Sterile Swabs x25, Test Tube with Buffer Solution x25, Test Card x25, Controls x2, Tube Holder x5, Instructions For Use x2 (1 English & 1 Spanish)
EUA Number	EUA230037