


INDICAID[®] OTC

UTI Test Rapid Test Strip


AVAILABLE IN
1-PACK AND 3-PACK



Urinary Tract Infection Prevalence



About 10 in 25 women and 3 in 25 men will have symptoms of a UTI during their lifetime.



In the United States, 40% of women will develop a UTI in their lifetime, with about 10% getting one each year.

8
Million

UTIs are the second most common type of infection in the body and are the reason for more than 8 million visits to the doctor each year.

Who is it for?

- Women are more prone to urinary tract infections due to their shorter urethra, which allows bacteria easier access to the bladder.
- Older adults are at higher risk of UTIs due to weakened immune systems, catheter use, and incomplete bladder emptying.
- Sexual activity can introduce bacteria into the urinary tract.

Test Type	Leukocytes (white blood cells) and Nitrite (OTC)
Sample Type	Urine
Time to Result	2 minutes
UTI Symptoms	Burning or pain during urination, frequent urination, fever, and back or groin pain.
Kit Contents	Test strips, instructions for use
HSA/FSA Approved	Yes
FDA 510(k)?	Yes

CONTACT US: 1-877-934-9344 OR USSALES@PHASESCI.COM

PHASESCIENTIFIC