

# INDICAID<sup>®</sup>

## Fentanyl Drug Detection Test

AVAILABLE IN  
25 TEST PACK



## Fentanyl Abuse Prevalence

# 74K

In 2023, the CDC estimated that fentanyl overdoses caused 74,702 deaths in the United States.



According to the National Institute on Drug Abuse over 115 million pills containing illicit fentanyl were seized by law enforcement in 2023.



Two milligrams of fentanyl can be lethal depending on a person's body size, tolerance, and past usage.

## Who is it for?

- Healthcare providers managing patients suspected of fentanyl use, to quickly confirm or rule out fentanyl exposure.
- Clinicians working in addiction recovery programs, to monitor and assess fentanyl presence in patients undergoing treatment.
- Law enforcement professionals conducting investigations or handling criminal justice cases involving suspected drug use, to provide rapid on-site confirmation of fentanyl presence.

<b>Test Type</b>	Rapid Lateral Fluid Immunoassay for the detection of Fentanyl in human urine
<b>Sample Type</b>	Urine Sample
<b>Time to Result</b>	5 minutes
<b>Symptoms of Fentanyl Overdose</b>	Decreased level of consciousness, slowed breathing, lack of response to stimulation, and constricted pupils.
<b>Kit Contents</b>	Test devices packaged individually in a foil pouch with desiccant and package insert.
<b>Urine Cutoff</b>	1 ng/mL
<b>FDA 510(k)?</b>	Yes