

INDICAID[®]

Multi-Drug Urine Test Cups

AVAILABLE IN
25 CUPS PER KIT



Drug Abuse Prevalence



In 2023, Over 107,543 overdose deaths were reported in the United States.

75%

Opioids are linked to about 75% of all overdose deaths in the United States.



About 22 million Americans have a substance use disorder, but only 10% receive treatment.

Who is it for?

- **Healthcare providers in primary and urgent care settings** use multi-drug test cups to quickly screen for various substances.
- **Clinicians working in addiction and rehabilitation programs** utilize these tests to monitor patients' substance use.
- **Medication management** utilize multi-drug test cups to monitor patients' adherence and detect potential drug misuse.

Test Type	Multi-drug urine test cups is a competitive binding, lateral flow immunochromatographic assays for the detection of multiple drugs of abuse (POC)
Substances Tested For	5 Panel: Amphetamines (AMP), Cocaine (COC), Methamphetamine (mAMP), Cannabis (THC), Opiates (OPI) 12 Panel: Amphetamines (AMP), Barbiturates (BAR), Buprenorphine (BUP), Benzodiazepines (BZO), Cocaine (COC), Methamphetamine (mAMP), MDMA/Ecstasy (MDMA), Morphine (MOP), Methadone (MTD), Oxycodone (OXY), Phencyclidine (PCP), Cannabis (THC)
Sample Type	Urine Sample
Time to Result	5 minutes
Kit Contents	25 multi-drug urine test cups in individual pouches, 1 package insert, 5 adulteration color comparison cards, 25 security seals, 25 gloves
Regulatory	FDA 510(k) & CLIA Waived